

POLYTECH Health & Aesthetics Bustiers

Postoperative Compression Therapy

POLYTECH

Health & Aesthetics



POLYTECH Health & Aesthetics Bustiers for postoperative compression therapy improve the aesthetic result of breast surgery.

POLYTECH Health & Aesthetics Bustiers

Postoperative Compression Therapy

Product Description

POLYTECH Health & Aesthetics Bustiers are available in 7 sizes with or without cups and in the colour black or white. They provide an optimal postoperative care after mammary augmentation, reconstruction or reduction.

The bustiers are made of mixed fibres with a high cotton content for greater patient comfort. They can be washed by hand at 30°C. Adjustable velcros, front opening with zip and three hooks and eyes guarantee an easy handling.

Advantages

- Minimizes the occurrence of hematoma
- Reduces postoperative oedema
- Reduces the occurrence of hypertrophic scars
- Holds implants and bandages in their position
- Facilitates wound inspection and handling

Important Information

POLYTECH Health & Aesthetics Bustiers are provided individually packed and non-sterile.

Attention: Please note the detailed instructions for use included with each product.



Order No.	Order No.	Underbust measurement cm	inch	Size
-----------	-----------	-----------------------------	------	------

POLYTECH Health & Aesthetics Bustiers, with cups

White	Black			
2110 P	2120 P	< 84	32	Petit
2110 S	2120 S	84 — 89	34	Small
2110 M	2120 M	89 — 94	36	Medium
2110 L	2120 L	94 — 104	38	Large
2110 XL	2120 XL	104 — 114	40	X-Large
2110 XXL	2120 XXL	114 — 124	42	XX-Large
2110 XXXL	2120 XXXL	124 — 134	44	XXX-Large

POLYTECH Health & Aesthetics Bustiers, without cups

White	Black			
2111 P	2121 P	< 84	32	Petit
2111 S	2121 S	84 — 89	34	Small
2111 M	2121 M	89 — 94	36	Medium
2111 L	2121 L	94 — 104	38	Large
2111 XL	2121 XL	104 — 114	40	X-Large
2111 XXL	2121 XXL	114 — 124	42	XX-Large
2111 XXXL	2121 XXXL	124 — 134	44	XXX-Large

Approximate measurements

POLYTECH Health & Aesthetics GmbH
 Altheimer Str. 32 • 64807 Dieburg • Germany
 Tel.: +49 (0) 60 71 98 63 0 • Fax: +49 (0) 60 71 98 63 30
 www.polytech-health-aesthetics.com • E-mail: info@polytechhealth.com

