



Malar implants



Nasal dorsal implants

Facial Implants

A selection of malar, nasal dorsal and chin implants for shaping an individual aesthetic facial contour



Chin implants

Implants made by POLYTECH – QUALITY made in Germany



Facial Implants, CE 0483

Information on Materials

Facial implants consist of a smooth silicone elastomer. The malar implants are equipped with pierced holes for fixation.

Production of the implants is based on a quality management according to ISO 9001 and ISO 13485. All materials used are of implantable medical quality.

Indications

- Augmentation
- Contour correction
- Reconstruction, e.g. after trauma
- Implant exchange

Important Information


Facial implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.

Order n°	Dimensions (cm)		
	A	B	C


Malar Implant ①			
10012-039	3.9	2.0	0.3
10012-042	4.2	2.2	0.34
10012-142	4.2	2.3	0.4

Malar Implant, anatomical design ②				
10011-050R	5.0	2.3	0.38	right
10011-050L	5.0	2.3	0.38	left
10011-052R	5.2	2.8	0.38	right
10011-052L	5.2	2.8	0.38	left
10011-055R	5.5	3.3	0.38	right
10011-055L	5.5	3.3	0.38	left
10011-056R	5.6	3.6	0.38	right
10011-056L	5,6	3,6	0,38	left




Order n°	Dimensions (cm)				
	A	B	C	D	E

Nasal Dorsal Implant					
10022-049	4.9	—	0.8	0.6	0.3
10022-149	4.9	—	0.8	0.54	0.5
10022-058	5.8	—	0.94	0.6	0.56




Nasal Dorsal Implant L-Shape						
10021-051	5.1	2.1	1.0	0.7	0.4	①
10021-052	5.2	2.1	0.84	0.65	0.4	②
10021-056	5.6	2.5	0.8	0.7	0.6	③




Order n°	Dimensions (cm)				
	A	B	C	D	E


Chin Implant, anatomical design 1					
10001-058	5.8	1.1	0.6	—	9.1
10001-059	5.9	1.2	0.6	—	10.1
10001-060	6.0	1.3	0.8	—	10.9
10001-062	6.2	1.3	0.8	—	11.4




Chin Implant, anatomical design 2					
10002-059	5.3	1.3	0.6	—	8.0
10002-073	7.3	1.4	0.6	—	9.4
10002-075	7.5	1.4	0.7	—	12.6



Chin Implant, anatomical design 3					
10003-057	5.7	1.3	0.8	—	8.2
10003-070	7.0	1.3	1.1	—	10.2



Chin Implant, anatomical design 4					
10004-154	5.4	1.5	0.3	—	6.1
10004-354	5.4	1.5	0.5	—	6.0
10004-554	5.4	1.5	0.7	—	6.1



approximate measurements

